

[implantswiss.com](http://implantswiss.com)

**NOVODENT SA**  
Y-Parc Technopole  
Avenue des Sciences 11  
CH-1400, Yverdon-les-Bains  
Switzerland

06.2020 Rev.03

 **implant**swiss  
SMART  
CHOICE

IMPLANT TECHNOLOGY

## » Content

Implantswiss Dental System	2
Swiss Precision	4
Product Traceability	5
Glass Package	6
Platform Color Coding	7
<b>» Implantswiss System Specifications</b>	8
S.R.A. Surface	12
Hybrid Design	13
Conical Connection with Hexagonal Index	14
Platform Switching	14
Self Cutting Edge	15
Hygienic Neck Design	15
<b>» Implantswiss Implants</b>	16
Bone Level Implant	18
Tissue Level Implant	19
<b>» Surgery</b>	20
Implantswiss Surgical Drills and Protocols	22
<b>» Short Implant</b>	28
Short Implant	29
Short Implant Surgical Kit	31
Short Implant Surgical Drills and Protocols	32
<b>» Zygoma Implant</b>	36
Zygoma Implant Surgical Kit	39
<b>» Implantswiss Prosthetic Components</b>	40
<b>» Bone Level System Prosthetic Components</b>	42
Bone Level Healing Abutments	45
Couple System	45
Solid System	47
Octa System	48
O-ring System	49
Locator System	50
Multi-Unit System	51
Bar System	52
<b>» Tissue Level System Prosthetic Components</b>	53
Tissue Level Healing Abutments	56
Couple System	56
Solid System	57
Octa System	58
O-ring System	59
Locator System	60
<b>» Zygoma Dental Implant System Prosthetic Components</b>	61
<b>» Digital Solutions</b>	62
Bone Level Digital Solutions	64
Tissue Level Digital Solutions	65
Multi Unit Abutment Digital Solutions	66
Octa Abutment Digital Solutions	66
Ti-Base Abutment	67
<b>» Index Code</b>	68

implantswiss

Switzerland

is at the heart of

**science, mechanical design  
and craftsmanship**



The well respected and acknowledged world authority...

Switzerland

has been one of the most important centers for craftsmanship that requires

**high precision and  
expertise** throughout history.

The products with impeccable mechanical designs, can be used for a lifetime and add value to your life.

**Implantswiss Dental System** owes its **excellent** structure to the inspiration, experience and value from history.

**Implantswiss** is a strong and reliable brand in the dental industry with its innovative solutions, unique surface design and high quality.

Implantswiss is devoted to Research and Development through cooperations with universities and crucial investments in clinical research. Implantswiss follows the latest trends and initiates extensive and contemporary projects that may lead the future of implant dentistry.

**Implantswiss  
Dental System**


is the most advanced stage  
in implant dentistry



# SWISS PRECISION



Implantswiss products are manufactured with the state-of-the-arts technologies and high precision. Implantswiss never compromises on its quality and reliability; the success and high quality of Implantswiss lies in this perfectionist approach in each step of the manufacturing process.




# PRODUCT TRACEABILITY



Each Implantswiss product can be easily traced from production to delivery.

Each product package has a barcode label with a lot number created by the ERP software system. All steps from production to stock, delivery and after sales can be monitored with the barcode system; therefore, the follow-up process can be conducted flawlessly.



# GLASS PACKAGE

## GLASS PACKAGE

The nature and human health have always been the primary concern for Implantswiss. Implantswiss, gives great importance to protection of nature and human health. Therefore, all Implantswiss Implants are packed in 100% natural and hygienic glass vials.

Since the glass structure doesn't contain any foreign material, glass packages are especially preferred for preserving the high value products.

The high barrier property, ease of sterilisation and unreactivity in chemical substances are among the many advantages of glass packages.

## SMART PACKAGE

The implant is simply removed from the glass bottle and placed in the osteotomy using the fixture driver.

Each product package has a label containing all necessary information such as the type, the length and the color code of the product.

The color coding and associated diameters of the implants are classified as given.



# PLATFORM COLOR CODING




## TISSUE LEVEL

## BONE LEVEL



» **IMPLANTSWISS  
SYSTEM  
SPECIFICATIONS**




# » SRA\* Surface


Following sandblasting, the surface of the implants is etched using Double Acid Etching technique.

Implantswiss S.R.A. surface has a homogenous surface roughness that perfectly mimics the bone morphology.

\* Sandblasted Rough Acid-etched Surface




### 3D-Area Parameters



Detail of the thread Bank and representation of the surface shape after removal of a 4th order polynomial

ISO 25178				
Height Parameters				
Sa	1,80	µm		
Sq	2,44	µm		
Feature Parameters				
S10z	15,1	µm	pruning = 5%	Ten point height
S5p	8,59	µm	pruning = 5%	Five point peak height
S5v	6,49	µm	pruning = 5%	Five point of height




Detail of the thread valley and representation of the surface shape after removal of a 4th order polynomial

ISO 25178				
Height Parameters				
Sa	1,25	µm		
Sq	1,60	µm		
Feature Parameters				
S10z	9,16	µm	pruning = 5%	Ten point height
S5p	3,39	µm	pruning = 5%	Five point peak height
S5v	5,76	µm	pruning = 5%	Five point of height

» Tapered form + Straight form =

# Hybrid Design

It is a unique implant design with cylindrical (straight) design on the neck area and conical (tapered) design through the apex of the implant.



The first 3 mm of Implantswiss implants running straight downward from the neck constitutes the straight form. The remaining portion constitutes the tapered form at a certain angle, depending on the length of the implant. The optimum combination of these two specifications creates the hybrid form. One advantage of the hybrid


design is to exert minimal stress on crestal cortical bone because of the straight portion at the neck, leading to an optimal load absorption rate. Due to its tapered form, 1/3 of the implant length is directly positioned in the osteotomy, which facilitate easy insertion of the implant in place.

## » Conical connection with hexagonal index

Creates an excellent implant-abutment connection

Implantswiss implants create an excellent implant-abutment interface, which prevents the formation of micro-gap and micro-leakage thanks to the universal conical connection and hexagonal index.

In Implantswiss Dental Implant System, all abutments are compatible with all fixtures with any diameter (except for Q 3.3) making a user friendly implant system for dentists.



## » Platform Switching

Protects the alveolar bone around the neck area

Micro-gap can be formed at the implant-abutment interface in Bone Level implants especially if they have a butt joint connection. Bacteria and by-products in these areas can cause serious problems such as bone resorption and periimplantitis.

Even though Implantswiss Implants have bacteria sealed conical connection, Switch Platform system further protects cervical bone by moving the implant-abutment interface away from alveolar bone. Implantswiss Implants with Switch Platform System are safer against bacteria and microleakage and have a longer life span.



## » Self Cutting Edge

Unique thread design that does not transfer stress to the bone

Implantswiss implants have a thread design which does not transfer stress to the bone during insertion. The cutting edges and self cutting feature provide ideal primary stability in any bone type without transferring stress thanks to its semi-aggressive thread design.




## » Hygienic Neck Design

Easy to maintain the hygiene

There are no micro threads on the neck area of Implantswiss implants therefore it is easier to maintain the hygiene when required.

Bacteria can easily colonize in the micro threads. In case of supra-alveolar positioning of the implant neck, the bone resorption risk is therefore minimized.





» IMPLANTSWISS  
IMPLANTS




# BONE LEVEL IMPLANT



Bone Level implants are designed to provide excellent outcome in any bone quality and in any location.

Implantswiss Bone Level implants allow to design an ideal smile in the aesthetic zone.




**System** » Bone Level  
**Material** » Ti Grade 4  
**Driver** » Bone Level driver,  
 Mini Bone Driver for 3.3 implant  
**Torque Value** » 35 Ncm

# TISSUE LEVEL IMPLANT



An excellent soft and hard tissue protection can be achieved with the smooth neck and location of implant-abutment interface of the Tissue Level implant.


The implant is designed to provide excellent outcomes in any case where the aesthetic outcome is not the primary concern.



**System** » Tissue Level  
**Material** » Ti Grade 4  
**Driver** » Tissue level driver  
**Torque Value** » 35 Ncm

## » SURGERY

### BONE LEVEL - TISSUE LEVEL




# IMPLANTSWISS FIXTURE DRILLS & PROTOCOLS




## SURGICAL DRILLS & PROTOCOLS

Self stoppered surgical drills eliminate the risk of complications both during the preparation of the osteotomy and insertion of the implant. Implantswiss drills are recommended to be used at 200-300 rpm to prevent over heating and reduce mechanical trauma, which will shorten the osseointegration time of the implant. The drills are designed to harvest autogenous bone when used below 50rpm without irrigation to be used as a graft material in regenerative procedures.

Implantswiss Surgical Kit contains a specific drill for each implant available in the system.




Implantswiss Surgical Kit includes various drills with different diameters and lengths. These drills are designed to simplify the surgical procedures resulting in shortened surgery time and decreased stress for the patients.




**ø3.7  
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø4.3  
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø4.8  
Fixture**




D1 Bone



D2-3 Bone




D4 Bone




**ø5.5  
Fixture**




D1 Bone



D2-3 Bone



D4 Bone



» **IMPLANTSWISS  
SHORT IMPLANTS**



# SHORT IMPLANT

In cases of insufficient bone volume, either above the mandibular canal or below the maxillary sinuses, short implants may be an alternative treatment to advanced surgical interventions aiming to regenerate hard tissues.

Due to aggressive thread design, Bone level short implants can be preferred mainly in maxilla, where the alveolar bone is relatively soft, whereas Tissue level short implant can be a good choice especially for the lower jaw.

## BONE LEVEL

	ø4.3	ø4.8	ø5.5
Length	ø4.3	ø4.8	ø5.5
5.0 mm	S-BMFSR4305	S-BMFSR4805	S-BMFSR5505
6.0 mm	S-BMFSR4306	S-BMFSR4806	S-BMFSR5506


## TISSUE LEVEL

	ø4.3	ø4.8	ø5.5
Length	ø4.3	ø4.8	ø5.5
4.0 mm	S-TFSR4304	S-TFSR4804	S-TWFSR5504
5.0 mm	S-TFSR4305	S-TFSR4805	S-TWFSR5505
6.0 mm	S-TFSR4306	S-TFSR4806	S-TWFSR5506

GH > 265 mm

# SHORT IMPLANT SURGICAL KIT


The surgical kit is recommended to be used with both Bone and Tissue Level short implants.






# SHORT IMPLANT DRILLS & PROTOCOLS

Implantswiss Short Implant Surgical Kit contains a specific drill for each implant available in the system.



Implantswiss Short Implant Surgical Kit includes various drills with different diameters and lengths. These drills are designed to simplify the surgical procedures resulting in shortened surgery time and decreased stress for the patients.




ø4.3  
Fixture




## BONE LEVEL

D4 / D3 Bone



## TISSUE LEVEL

D2 / D1 Bone




Bone Level Short Implants are recommended to be used only in D3/D4 and Tissue Level Implants only in D1/D2 bone.

**ø4.8  
Fixture**



BONE LEVEL

D4 / D3 Bone




**ø5.5  
Fixture**




BONE LEVEL

D4 / D3 Bone




TISSUE LEVEL

D2 / D1 Bone



TISSUE LEVEL

D2 / D1 Bone



» **IMPLANTSWISS  
ZYGOMA IMPLANTS**



# ZYGOMA IMPLANT

Zygomatic implants are designed for cases where traditional implants can not be used and for patients who want to benefit from immediate loading protocols even in the presence of minimum residual maxillary bone height.

Implantswiss Zygomatic Implant offers three different product options.

## ZYGOMA IMPLANT

ø4.3



35 mm	S-ZYG35
37.5 mm	S-ZYG37.5
40 mm	S-ZYG40
42.5 mm	S-ZYG42.5
45 mm	S-ZYG45
47.5 mm	S-ZYG47.5
50 mm	S-ZYG50
52.5 mm	S-ZYG52.5
55 mm	S-ZYG55
57.5 mm	S-ZYG57.5
60 mm	S-ZYG60

## ATROPHIC MAXILLARY ZYGOMA IMPLANT

ø4.3



35 mm	S-ZAM35
37.5 mm	S-ZAM37.5
40 mm	S-ZAM40
42.5 mm	S-ZAM42.5
45 mm	S-ZAM45
47.5 mm	S-ZAM47.5
50 mm	S-ZAM50
52.5 mm	S-ZAM52.5
55 mm	S-ZAM55
57.5 mm	S-ZAM57.5
60 mm	S-ZAM60


## ONCOLOGICAL ZYGOMA IMPLANT

ø4.3



25 mm	SZ-ONC25
27.5 mm	SZ-ONC27.5
30 mm	SZ-ONC30
32.5 mm	SZ-ONC32.5

# ZYGOMA IMPLANT SURGICAL KIT



## SURGICAL DRILLS

Point drill    Pilot drill 3.5    Pilot drill Long    Twist drill 2.5 Short    Twist drill 2.5 Long    Twist drill 2.9 Short    Twist drill 2.9 Long    Twist drill 3.5 Short    Twist drill 3.5 Long

















» **IMPLANTSWISS  
PROSTHETIC COMPONENTS**




















# BONE LEVEL















## COUPLE SYSTEM

FLOW	FIXTURE	HEALING	IMPRESSION COPING	LAB ANALOG	ABUTMENT	SCREW
		 Closing Screw	 Pick-up Hex		 Couple Hex	
		 Healing Abutment	 Transfer Hex		 Angled 15° 25°	
TOOL	 Ratchet Driver	 Hex Screw Driver	 Hex Screw Driver			 Hex Screw Driver













## OCTA SYSTEM

FLOW	FIXTURE	HEALING	SUB OCTA ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER	SCREW
		 Closing Screw	 Octa Healing Cap	 Pick-up		 Octa	
		 Healing Abutment	 Octa Abutment	 Transfer		 N-Octa	
TOOL	 Ratchet Driver	 Hex Screw Driver	 Octa Driver			 Octa/N-Octa	 Hex Screw Driver

## SOLID SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION CAP	LAB ANALOG	PLASTIC COPING
		 Closing Screw	 Solid			 Single
		 Healing Abutment	 Cap			 Bridge
TOOL	 Ratchet Driver	 Hex Screw Driver	 Solid Driver	 Hex Screw Driver		 Hex Screw Driver

## O-RING SYSTEM

FLOW	FIXTURE	HEALING	O-RING ABUTMENT	O-RING LAB ANALOG	RETAINER
		 Closing Screw			 S-ORR
		 Healing Abutment			 S-ORRc
TOOL	 Ratchet Driver	 Hex Screw Driver	 O-Ring Driver		 S-ORRo

## LOCATOR SYSTEM

FLOW	FIXTURE	HEALING	LOCATOR ABUTMENT	LOCATOR LAB ANALOG	INSERTER
TOOL	 Ratchet Driver	 Hex Screw Driver	 Locator Driver		 S-LRI

## MULTI-UNIT SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER
TOOL	 Ratchet Driver	 Hex Screw Driver	 Hex Screw Driver Multi-Unit Driver	 Hex Screw Driver		 Hex Screw Driver

## BONE LEVEL HEALING ABUTMENTS

### BONE LEVEL HEALING ABUTMENTS

	 ø3.7 <small>only for ø3.3 implant</small>	 ø4.5	 ø5.5	 ø6.5
1 mm	-	S-BHA4510	S-BHA5510	S-BHA6510
2 mm	-	S-BHA4520	S-BHA5520	S-BHA6520
2.5 mm	S-BMHA3725	-	-	-
3 mm	-	S-BHA4530	S-BHA5530	S-BHA6530
3.5 mm	S-BMHA3735	-	-	-
4 mm	-	S-BHA4540	S-BHA5540	S-BHA6540
5 mm	S-BMHA3750	S-BHA4550	S-BHA5550	S-BHA6550
7 mm	S-BMHA3770	-	-	-

System » Bone Level  
Material » Ti Grade 5 ELI  
Driver » Hex Screw Driver

### BONE LEVEL SHORT IMPLANT HEALING ABUTMENTS

	 ø5.0	 ø6.0	 ø8.0
2 mm	S-BMHA5020	S-BMHA6020	S-BMHA8020
3 mm	S-BMHA5030	S-BMHA6030	S-BMHA8030

System » Bone Level  
Material » Ti Grade 5 ELI  
Driver » Hex Screw Driver




## BONE LEVEL COUPLE SYSTEM

### BONE LEVEL COUPLE ABUTMENTS

	 ø3.7 <small>only for ø3.3 implant</small>	 ø4.5	 ø5.5	 ø6.5
1 mm	S-BMCA37107	S-BCA45107	S-BCA55107	S-BCA65107
2 mm	S-BMCA37207	S-BCA45207	S-BCA55207	S-BCA65207
3 mm	S-BMCA37307	S-BCA45307	S-BCA55307	S-BCA65307
4 mm	S-BMCA37407	S-BCA45407	S-BCA55407	S-BCA65407
5 mm	S-BMCA37507	S-BCA45507	S-BCA55507	S-BCA65507





System » Bone Level  
Material » Ti Grade 5 ELI  
Torque Value » 25 Ncm  
Driver » Hex Screw Driver

## BONE LEVEL SHORT IMPLANT COUPLE ABUTMENTS

			
	Ø5.0	Ø6.0	Ø8.0
<b>2 mm</b>	S-BMCA5020	S-BMCA6020	S-BMCA8020
<b>3 mm</b>	S-BMCA5030	S-BMCA6030	S-BMCA8030




**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver

## BONE LEVEL ANGLED ABUTMENT

				
	Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5 <small>15° 25°</small>	Ø5.5 <small>15° 25°</small>	Ø6.5 <small>15° 25°</small>
<b>1 mm (15°)</b>	S-BMAA37151	<b>2 mm 15°</b> S-BAA45152	S-BAA55152	S-BAA65152
<b>2 mm (15°)</b>	S-BMAA37152	<b>2 mm 25°</b> S-BAA45252	S-BAA55252	S-BAA65252
<b>4 mm (15°)</b>	S-BMAA37154	<b>4 mm 15°</b> S-BAA45154	S-BAA55154	S-BAA65154
		<b>4 mm 25°</b> S-BAA45254	S-BAA55254	S-BAA65254





**System** » Bone Level // **Material** » Ti Grade 5 ELI // **Torque Value** » 25 Ncm // **Driver** » Hex Screw Driver

## TEMPORARY ABUTMENTS

		
Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5	Ø5.5
S-BMTA37	S-BTA45	S-BTA55



**System** » Bone Level  
**Material** » PEEK  
**Torque Value** » 20 Ncm  
**Driver** » Hex Screw Driver

## BONE LEVEL IMPRESSION COPING PICK UP

			
Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5	Ø5.5	Ø6.5
S-BMICP37L S-BMICP37S	S-BICP45L S-BICP45S	S-BICP55L S-BICP55S	S-BICP65L S-BICP65S



**System** » Bone Level  
**Material** » Ti Grade 5 ELI

## BONE LEVEL SHORT IMPLANT IMPRESSION COPING PICK UP

	
Ø5.0	Ø6.0
S-BMSICP50	S-BMSICP60

**System** » Bone Level  
**Material** » Ti Grade 5 ELI

## BONE LEVEL IMPRESSION COPING TRANSFER

	
Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5
S-BMICT37L S-BMICT37S	S-BIP45



**System** » Bone Level  
**Material** » Ti Grade 5 ELI

## BONE LEVEL SHORT IMPLANT IMPRESSION COPING TRANSFER

	
S-BMSICT50	S-BMSICT60

**System** » Bone Level  
**Material** » Ti Grade 5 ELI




## BONE LEVEL FIXTURE LAB ANALOG

	
Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5
S-BMFLA	S-BFLA

**System** » Bone Level  
**Material** » Ti Grade 5 ELI

# BONE LEVEL SOLID SYSTEM


## BONE LEVEL SOLID ABUTMENTS

			
	Ø3.7 <small>only for ø3.3 implant</small>	Ø4.5	Ø5.5
<b>1 mm</b>	S-BMSA3710	S-BSA4510	S-BSA5510
<b>2 mm</b>	S-BMSA3720	S-BSA4520	S-BSA5520
<b>3 mm</b>	S-BMSA3730	S-BSA4530	S-BSA5530
<b>4 mm</b>	S-BMSA3740	S-BSA4540	S-BSA5540
<b>5 mm</b>	S-BMSA3750	S-BSA4550	S-BSA5550

**System** » Bone Level  
**Material** » Ti Grade 5 ELI  
**Torque Value** » 25 Ncm  
**Driver** » Hex Screw Driver




## BONE LEVEL SOLID ABUTMENT CAP




**System »** Bone Level Solid System  
**Material »** POM

## BONE LEVEL SOLID ABUTMENT IMPRESSION CAP



**System »** Bone Level Solid System  
**Material »** Plastic

## BONE LEVEL SOLID LAB ANALOG



**System »** Bone Level Solid System  
**Material »** Ti Grade 5 ELI


# BONE LEVEL OCTA SYSTEM

## BONE LEVEL OCTA ABUTMENTS



**System »** Bone Level  
**Material »** Ti Grade 5 ELI  
**Torque Value »** : Abutment Screw  
 25 Ncm Occlusal Screw 20 Ncm  
**Driver »** Octa Driver

## BONE LEVEL OCTA ABUTMENT HEALING CAP



**System »** Octa System  
**Material »** Ti Grade 4

## BONE LEVEL OCTA SYSTEM PLASTIC CYLINDER



**System »** Octa System  
**Material »** POM

## BONE LEVEL OCTA SYSTEM TITANIUM CYLINDER



**System »** Octa System  
**Material »** Ti Grade 5 ELI

## BONE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP



## BONE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER



**System »** Octa System  
**Material »** Ti Grade 5 ELI


## BONE LEVEL OCTA ABUTMENT LAB ANALOG



**System »** Octa System  
**Material »** Ti Grade 5 ELI

# BONE LEVEL O-RING SYSTEM

## BONE LEVEL O-RING ABUTMENTS



**System »** Bone Level  
**Material »** Ti Grade 5 ELI  
**Torque Value »** 25 Ncm  
**Driver »** O-Ring Driver

## DIRECTIONAL RINGS



S-ORRa

## O-RING LAB ANALOG



S-ORLA

System » O-Ring System  
Material » Ti Grade 5 ELI

## RETAINER



		PINK Soft 900 Gr	YELLOW ExtraSoft 500 Gr	GREEN Elastic 350 Gr
S-ORR	S-ORRc	S-ORG-Pcap	S-ORG-Ycap	S-ORG-Gcap

## BONE LEVEL LOCATOR SYSTEM

### BONE LEVEL MINI LOCATOR ABUTMENTS

only for ø3.3 implant

1.5 mm	S-BMLA15
2.5 mm	S-BMLA25
3.5 mm	S-BMLA35
4.5 mm	S-BMLA45
5.5 mm	S-BMLA55
6.5 mm	S-BMLA65

### BONE LEVEL LOCATOR ABUTMENTS

1.5 mm	S-BLA15
2.5 mm	S-BLA25
3.5 mm	S-BLA35
4.5 mm	S-BLA45
5.5 mm	S-BLA55
6.5 mm	S-BLA65



System » Bone Level  
Material » Ti Grade 5 ELI  
Torque Value » 25 Ncm  
Driver » Locator Driver

## BONE LEVEL LOCATOR LAB ANALOG




S-BLLA

## LOCATOR MATRIX



TITANIUM HOUSING CAP


S-LMht




BLACK CAP  
Processing




VIOLET CAP  
Rigid Retention 2.7 Kg



PINK CAP  
Soft Retention 1.2 Kg



CLEAR CAP  
Standard Retention 1.8 Kg



YELLOW CAP  
Extra Soft Retention 0.6 Kg

S-LRI

## BONE LEVEL MULTI UNIT SYSTEM

### MULTI UNIT ABUTMENTS



ø4.8

only for ø3.3 implant

0.5 mm	S-BMMUA4805
1.5 mm	S-BMMUA4815
2.5 mm	S-BMMUA4825
3.5 mm	S-BMMUA4835
4.5 mm	S-BMMUA4845



ø4.8

0.5 mm	S-BMUA4805
1.5 mm	S-BMUA4815
2.5 mm	S-BMUA4825
3.5 mm	S-BMUA4835
4.5 mm	S-BMUA4845

System » Bone Level  
Material » Ti Grade 5 ELI  
Torque Value » : Abutment Screw  
25 Ncm Oclusal Screw 20 Ncm  
Driver » Multi Unit Driver

### MULTI UNIT ANGLED ABUTMENTS



ø4.8 (17°)

1.0 mm	S-BMUA48171
2.5 mm	S-BMUA48172
3.5 mm	S-BMUA48173



ø4.8 (35°)

1.0 mm	S-BMUA48351
2.5 mm	S-BMUA48352
3.5 mm	S-BMUA48353

System » Bone Level  
Material » Ti Grade 5 ELI  
Torque Value » : Abutment Screw  
25 Ncm Oclusal Screw 20 Ncm  
Driver » 1,2 Hex Driver

### MULTI UNIT ABUTMENT HEALING CAP



ø4.8

4 mm	S-MUAHCap
6 mm	S-MUAHCapL

System » Multi Unit System  
Material » Ti Grade 5 ELI

### MULTI UNIT IMPRESSION COPING PICK UP



ø4.8

S-BMUAICP

System » Multi Unit System  
Material » Ti Grade 5 ELI

## MULTI UNIT ABUTMENT PLASTIC CYLINDER



S-MUAPICyN

**System** » Multi Unit System  
**Material** » POM

## MULTI UNIT ABUTMENT TITANIUM CYLINDER



S-MUATICyN

**System** » Multi Unit System  
**Material** » Ti Grade 5 ELI

## MULTI UNIT LAB ANALOG



S-MUALA  
Ø 4.8×14mm

**System** » Bone Level  
**Material** » Ti Grade 5 ELI

## BAR SYSTEM

### BAR BURNOUT

S-IB



### BAR CLIP YELLOW

S-IBClip - y  
Medium Retention 1.8 Kg



### TOOL FOR INSERTING CLIP



S-BMIA

### BAR POSITIONING

S-IBP



### BAR CLIP PINK

S-IBClip - p  
Soft Retention 1.2 Kg



### KEY FOR PARALLELOMETER



S-TFIC

### BAR HOUSING

S-IBH



### BAR CLIP RED

S-IBClip - r  
Extra Soft Retention 0.8 Kg



## TISSUE LEVEL
















### COUPLE SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	LAB ANALOG	ABUTMENT	SCREW
 ø3.7 ø4.3-L06 ø4.3  ø4.8-L06 ø4.8 ø5.5	 Cover Screw	 Pick-up Octa	 Transfer Octa		 Couple Octa  Angled 15° 25°  Synocta Octa	
TOOL	 Ratchet Driver	 Hex Screw Driver	 Hex Screw Driver			 Hex Screw Driver






### SOLID SYSTEM

FLOW	FIXTURE	HEALING	ABUTMENT	IMPRESSION CAP	LAB ANALOG	PLASTIC COPING
 ø3.7 ø4.3-L06 ø4.3  ø4.8-L06 ø4.8 ø5.5	 Closing Screw	 Healing Abutment	 Cap		 Lab Analog	 Crown  Bridge
TOOL	 Ratchet Driver	 Hex Screw Driver	 Solid Driver	 Hex Screw Driver		 Hex Screw Driver













## OCTA SYSTEM

FLOW	FIXTURE	HEALING	OCTA ABUTMENT	IMPRESSION COPING	LAB ANALOG	LAB CYLINDER	SCREW
 ø3.7 ø4.3-L06 ø4.3  ø4.8-L06 ø4.8 ø5.5		 Cover Screw  Healing Abutment		 Pick-up  Transfer		 Octa  N-Octa	
TOOL	 Ratchet Driver	 Hex Screw Driver	 Octa Driver			 Octa/N-Octa	 Hex Screw Driver

## LOCATOR SYSTEM

FLOW	FIXTURE	HEALING	LOCATOR ABUTMENT	LOCATOR LAB ANALOG	INSERTER
 ø3.7 ø4.3-L06 ø4.3  ø4.8-L06 ø4.8 ø5.5		 Closing Screw  Healing Abutment			 S-LMHt
TOOL	 Ratchet Driver	 Hex Screw Driver	 Locator Driver		 S-LRI

## O-RING SYSTEM

FLOW	FIXTURE	HEALING	O-RING ABUTMENT	O-RING LAB ANALOG	RETAINER
 ø3.7 ø4.3-L06 ø4.3  ø4.8-L06 ø4.8 ø5.5		 Closing Screw  Healing Abutment			 S-ORR  S-ORRc
TOOL	 Ratchet Driver	 Hex Screw Driver	 O-Ring Driver		 S-ORR

# TISSUE LEVEL HEALING ABUTMENTS

## TISSUE LEVEL HEALING ABUTMENTS



3.5 mm S-THA35 S-THAW35  
5.5 mm S-THA55 S-THAW55

only for ø5.5 implant

System » Tissue Level  
Material » Ti Grade 4  
Torque Value » 20 Ncm  
Driver » Hex Screw Driver

# TISSUE LEVEL COUPLE SYSTEM

## TISSUE LEVEL COUPLE ABUTMENTS



1 mm S-TCA10 S-TCAW10  
2 mm S-TCA20 S-TCAW20  
3 mm S-TCA30 S-TCAW30  
4 mm S-TCA40 S-TCAW40

only for ø5.5 implant

System » Tissue Level  
Material » Ti Grade 4  
Torque Value » 25 Ncm  
Driver » Hex Screw Driver

## TISSUE LEVEL ANGLED ABUTMENT



15° S-TAA15  
25° S-TAA25

System » Tissue Level  
Material » Ti Grade 4  
Torque Value » 25 Ncm  
Driver » Hex Screw Driver


## TISSUE LEVEL SYNOCTA ABUTMENT



6 mm S-TSOA

System » Tissue Level  
Material » Ti Grade 4  
Torque Value » 25 Ncm  
Driver » Hex Screw Driver

## TISSUE LEVEL IMPRESSION COPING PICK UP



OCTA S-TICP S-TICPW

only for ø5.5 implant

System » Tissue Level  
Material » Ti Grade 5 ELI

## TISSUE LEVEL IMPRESSION COPING TRANSFER



OCTA S-TICT

System » Tissue Level  
Material » Ti Grade 5 ELI

## TISSUE LEVEL FIXTURE LAB ANALOG



S-TFLA S-TFLLA

only for ø5.5 implant

System » Tissue Level  
Material » Ti Grade 5 ELI

# TISSUE LEVEL SOLID SYSTEM

## TISSUE LEVEL SOLID ABUTMENTS



ø3.5 S-TSA40 S-TSAL40  
ø4.5 S-TSAL50  
only for ø5.5 implant  
4 mm S-TSA40 S-TSAL40  
5 mm - S-TSAL50  
5.5 mm S-TSA55 -  
7 mm S-TSA70 S-TSAL70

System » Tissue Level  
Material » Ti Grade 4  
Torque Value » 25 Ncm  
Driver » Solid Driver

## TISSUE LEVEL SOLID ABUTMENT POSITIONING CYLINDER & IMPRESSION CAP



S-TSAPoCy S-TSAICap

System » Tissue Level Solid System  
Material » Plastic



## TISSUE LEVEL SOLID PLASTIC COPING



CROWN S-TSAPICc  
BRIDGE S-TSAPICb

System » Tissue Level Solid System  
Material » POM

## TISSUE LEVEL SOLID ABUTMENT LAB ANALOG

		
	only for ø5.5 implant	
<b>4 mm</b>	S-TSALA40	S-TSALLA40
<b>5.5 mm</b>	S-TSALA55	S-TSALLA55
<b>7 mm</b>	S-TSALA70	S-TSALLA70

**System** » Tissue Level Solid System  
**Material** » Ti Grade 5 ELI

## TISSUE LEVEL OCTA SYSTEM

### TISSUE LEVEL OCTA ABUTMENT



**System** » Tissue Level  
**Material** » Ti Grade 4  
**Torque Value** » : Abutment Screw 25 Ncm Oclusal Screw 20 Ncm  
**Driver** » Octa Driver

### TISSUE LEVEL OCTA ABUTMENT HEALING CAP

	
<b>4 mm</b>	S-OAHCap
<b>6 mm</b>	S-OAHCapL

**System** » Octa System  
**Material** » Ti Grade 5 ELI

### TISSUE LEVEL OCTA ABUTMENT PLASTIC CYLINDER OCTA / NON-OCTA

		
<b>OCTA</b>	S-OAPICy	
<b>NON-OCTA</b>	S-OAPICyN	

**System** » Octa System  
**Material** » POM

### TISSUE LEVEL OCTA ABUTMENT TITANIUM CYLINDER OCTA / NON-OCTA

	
<b>OCTA</b>	S-OATICy
<b>NON-OCTA</b>	S-OATICyN

**System** » Octa System  
**Material** » Ti Grade 5 ELI

## TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP OCTA / NON-OCTA

	
	ø5.5
<b>OCTA</b>	S-ICOAP
<b>N-OCTA</b>	S-ICOAPN

**System** » Octa System  
**Material** » Ti Grade 5 ELI

## TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER OCTA / NON-OCTA

	
	ø5.0
<b>OCTA</b>	S-ICOAT
<b>N-OCTA</b>	S-ICOATN

**System** » Octa System  
**Material** » Ti Grade 5 ELI


## TISSUE LEVEL OCTA LAB ANALOG



**System** » Octa System  
**Material** » Ti Grade 5 ELI

## TISSUE LEVEL O-RING SYSTEM

### TISSUE LEVEL O-RING ABUTMENTS

	
<b>0.5 mm</b>	S-TORA05
<b>2 mm</b>	S-TORA20
<b>4 mm</b>	S-TORA40

**System** » Tissue Level  
**Material** » Ti Grade 4  
**Torque Value** » 25 Ncm  
**Driver** » O-Ring Driver

### DIRECTIONAL RINGS



### O-RING LAB ANALOG



**System** » O-Ring System  
**Material** » Ti Grade 5 ELI

## RETAINER



PINK  
Soft 900 Gr

YELLOW  
ExtraSoft 500 Gr

GREEN  
Elastic 350 Gr

S-ORR

S-ORRc

S-ORG-Pcap

S-ORG-Ycap

S-ORG-Gcap

## TISSUE LEVEL LOCATOR SYSTEM

### TISSUE LEVEL LOCATOR ABUTMENTS



1.5 mm	S-TLA15
2.5 mm	S-TLA25
3.5 mm	S-TLA35
4.5 mm	S-TLA45
5.5 mm	S-TLA55
6.5 mm	S-TLA65

System » Tissue Level  
Material » Ti Grade 5 ELI  
Torque Value » 25 Ncm  
Driver » Locator Driver

### TISSUE LEVEL LOCATOR LAB ANALOG



S-TLLA

### LOCATOR MATRIX



S-LMHT

S-LRI

## ZYGOMA DENTAL IMPLANT SYSTEM

### ZYGOMA COUPLE ABUTMENT



3 mm	S-ZCA4330
5 mm	S-ZCA4350

System » Zygoma  
Material » Ti Grade 5 ELI  
Torque Value » 25 Ncm  
Driver » Hex Screw Driver

### ANGLED ABUTMENT



3 mm	S-ZAA43203
5 mm	S-ZAA43205

System » Zygoma  
Material » Ti Grade 5 ELI  
Torque Value » 25 Ncm  
Driver » Hex Screw Driver

### MULTI UNIT ABUTMENT



∅4.3 (17°)

0.5 mm	S-ZMUA431705
1.0 mm	S-ZMUA43171
2.0 mm	S-ZMUA43172
3.0 mm	S-ZMUA43173
5.0 mm	S-ZMUA43175

System » Zygoma  
Material » Ti Grade 5 ELI  
Torque Value » 25 Ncm  
Driver » Hex Screw Driver

### LAB ANALOG



S-ZLA

System » Zygoma // Material » Ti Grade 5 ELI

### HEALING CAP



3 mm	S-ZHA4330
5 mm	S-ZHA4350

System » Zygoma // Torque Value » 25 Ncm

### IMPRESSION COPING PICK UP



S-ZICP

System » Zygoma // Material » Ti Grade 5 ELI

### PLASTIC CYLINDER



S-ZPICYN

System » Zygoma // Material » POM

» **DIGITAL  
SOLUTIONS**





# BONE LEVEL

## SCAN BODY



only for ø3.3 implant  
S-PMSB



S-PSB

System » Digital System  
Driver » Hex Screw Driver

## DIGITAL ANALOG



only for ø3.3 implant  
S-BMFDLA



S-BFDLA

System » Digital System  
Driver » Hex Screw Driver

## PREMILL



only for ø3.3 implant  
ø10 mm S-BMPA10  
ø14 mm S-BMPA14  
ø16 mm -



S-BPA10  
S-BPA14  
S-BPA16

System » Digital System  
Material » Ti Grade 5 ELI  
Driver » Hex Screw Driver

# TISSUE LEVEL

## SCAN BODY



S-TSB



only for ø5.5 implant  
S-TWSB

System » Digital System  
Driver » Hex Screw Driver

## DIGITAL ANALOG



S-TFDLA



only for ø5.5 implant  
S-TWFDLA

System » Digital System  
Driver » Hex Screw Driver

## PREMILL



ø10 mm S-TPA10  
ø14 mm S-TPA14  
ø16 mm S-TPA16



only for ø5.5 implant  
S-TWPA10  
S-TWPA14  
S-TWPA16

System » Digital System  
Material » Ti Grade 5 ELI  
Driver » Hex Screw Driver

# MULTI UNIT ABUTMENT

## SCAN BODY



S-PMUSB

System » Digital System  
Driver » Hex Screw Driver

## DIGITAL ANALOG



S-BMDLA

System » Digital System  
Driver » Hex Screw Driver

## ZR BASE



S-MUTB

System » Digital System  
Driver » Hex Screw Driver

# OCTA ABUTMENT

## SCAN BODY



S-POSB

System » Digital System  
Driver » Hex Screw Driver

## DIGITAL ANALOG



S-BODLA

System » Digital System  
Driver » Hex Screw Driver

## ZR BASE



S-OTB

System » Digital System  
Driver » Hex Screw Driver

# Ti-BASE ABUTMENT

## Ti-BASE SCAN BODY



only for ø3.3 implant  
S-TIMSB



S-TISB

System » Digital System  
Driver » Hex Screw Driver

## Ti-BASE ABUTMENT



only for ø3.3 implant  
S-BMTIA05



S-BTIA05

System » Digital System  
Driver » Hex Screw Driver

# INDEX CODE

PRODUCT	PLATFORM	HEIGHT	CODE
---------	----------	--------	------

## BONE LEVEL FIXTURES

### BONE LEVEL MINI FIXTURE

Bone Level Mini Fixture Straight SRA	ø 3.3	8 mm	S-BMFSR3308
Bone Level Mini Fixture Straight SRA	ø 3.3	10 mm	S-BMFSR3310
Bone Level Mini Fixture Straight SRA	ø 3.3	12 mm	S-BMFSR3312
Bone Level Mini Fixture Straight SRA	ø 3.3	14 mm	S-BMFSR3314

### BONE LEVEL SHORT FIXTURE

Bone Level Short Fixture Straight SRA	ø 4.3	5 mm	S-BMFSR4305
Bone Level Short Fixture Straight SRA	ø 4.3	6 mm	S-BMFSR4306

Bone Level Short Fixture Straight SRA	ø 4.8	5 mm	S-BMFSR4805
Bone Level Short Fixture Straight SRA	ø 4.8	6 mm	S-BMFSR4806

Bone Level Short Fixture Straight SRA	ø 5.5	5 mm	S-BMFSR5505
Bone Level Short Fixture Straight SRA	ø 5.5	6 mm	S-BMFSR5506

### BONE LEVEL FIXTURE

Bone Level Wide Fixture Hybrid SRA	ø 3.7	8 mm	S-BFHR3708
Bone Level Wide Fixture Hybrid SRA	ø 3.7	10 mm	S-BFHR3710
Bone Level Wide Fixture Hybrid SRA	ø 3.7	12 mm	S-BFHR3712
Bone Level Wide Fixture Hybrid SRA	ø 3.7	14 mm	S-BFHR3714

Bone Level Wide Fixture Hybrid SRA	ø 4.3	8 mm	S-BFHR4308
Bone Level Wide Fixture Hybrid SRA	ø 4.3	10 mm	S-BFHR4310
Bone Level Wide Fixture Hybrid SRA	ø 4.3	12 mm	S-BFHR4312
Bone Level Wide Fixture Hybrid SRA	ø 4.3	14 mm	S-BFHR4314

Bone Level Wide Fixture Hybrid SRA	ø 4.8	8 mm	S-BFHR4808
Bone Level Wide Fixture Hybrid SRA	ø 4.8	10 mm	S-BFHR4810
Bone Level Wide Fixture Hybrid SRA	ø 4.8	12 mm	S-BFHR4812
Bone Level Wide Fixture Hybrid SRA	ø 4.8	14 mm	S-BFHR4814

### BONE LEVEL WIDE FIXTURE

Bone Level Wide Fixture Hybrid SRA	ø 5.5	8 mm	S-BWFHR5508
Bone Level Wide Fixture Hybrid SRA	ø 5.5	10 mm	S-BWFHR5510
Bone Level Wide Fixture Hybrid SRA	ø 5.5	12 mm	S-BWFHR5512

### HEALING ABUTMENTS

Healing Abutment	ø 3.7	2.5 mm	S-BMHA3725
Healing Abutment	ø 3.7	3.5 mm	S-BMHA3735
Healing Abutment	ø 3.7	5 mm	S-BMHA3750
Healing Abutment	ø 3.7	7 mm	S-BMHA3770

Healing Abutment	ø 4.5	1 mm	S-BHA4510
Healing Abutment	ø 4.5	2 mm	S-BHA4520

PRODUCT	PLATFORM	HEIGHT	CODE
---------	----------	--------	------

Healing Abutment	ø 4.5	3 mm	S-BHA4530
Healing Abutment	ø 4.5	4 mm	S-BHA4540
Healing Abutment	ø 4.5	5 mm	S-BHA4550

Healing Abutment	ø 5.5	1 mm	S-BHA5510
Healing Abutment	ø 5.5	2 mm	S-BHA5520
Healing Abutment	ø 5.5	3 mm	S-BHA5530
Healing Abutment	ø 5.5	4 mm	S-BHA5540
Healing Abutment	ø 5.5	5 mm	S-BHA5550

Healing Abutment	ø 6.5	1 mm	S-BHA6510
Healing Abutment	ø 6.5	2 mm	S-BHA6520
Healing Abutment	ø 6.5	3 mm	S-BHA6530
Healing Abutment	ø 6.5	4 mm	S-BHA6540
Healing Abutment	ø 6.5	5 mm	S-BHA6550

Mini Healing Abutment for Short Implant	ø 5.0	2 mm	S-BMHA5020
Mini Healing Abutment for Short Implant	ø 5.0	3 mm	S-BMHA5030

Mini Healing Abutment for Short Implant	ø 6.0	2 mm	S-BMHA6020
Mini Healing Abutment for Short Implant	ø 6.0	3 mm	S-BMHA6030

Mini Healing Abutment for Short Implant	ø 8.0	2 mm	S-BMHA8020
Mini Healing Abutment for Short Implant	ø 8.0	3 mm	S-BMHA8030

## BONE LEVEL COUPLE SYSTEM

### COUPLE ABUTMENTS

Mini Couple Abutment H7	ø 3.7	1 mm	S-BMCA37107
Mini Couple Abutment H7	ø 3.7	2 mm	S-BMCA37207
Mini Couple Abutment H7	ø 3.7	3 mm	S-BMCA37307
Mini Couple Abutment H7	ø 3.7	4 mm	S-BMCA37407
Mini Couple Abutment H7	ø 3.7	5 mm	S-BMCA37507

Couple Abutment H7	ø 4.5	1 mm	S-BCA45107
Couple Abutment H7	ø 4.5	2 mm	S-BCA45207
Couple Abutment H7	ø 4.5	3 mm	S-BCA45307
Couple Abutment H7	ø 4.5	4 mm	S-BCA45407
Couple Abutment H7	ø 4.5	5 mm	S-BCA45507

Couple Abutment H7	ø 5.5	1 mm	S-BCA55107
Couple Abutment H7	ø 5.5	2 mm	S-BCA55207
Couple Abutment H7	ø 5.5	3 mm	S-BCA55307
Couple Abutment H7	ø 5.5	4 mm	S-BCA55407
Couple Abutment H7	ø 5.5	5 mm	S-BCA55507

Couple Abutment H7	ø 6.5	1 mm	S-BCA65107
Couple Abutment H7	ø 6.5	2 mm	S-BCA65207
Couple Abutment H7	ø 6.5	3 mm	S-BCA65307
Couple Abutment H7	ø 6.5	4 mm	S-BCA65407
Couple Abutment H7	ø 6.5	5 mm	S-BCA65507

PRODUCT	PLATFORM	HEIGHT	CODE
Mini Couple Abutment for Short Implant	ø 5.0	2 mm	S-BMCA5020
Mini Couple Abutment for Short Implant	ø 5.0	3 mm	S-BMCA5030
Mini Couple Abutment for Short Implant	ø 6.0	2 mm	S-BMCA6020
Mini Couple Abutment for Short Implant	ø 6.0	3 mm	S-BMCA6030
Mini Couple Abutment for Short Implant	ø 8.0	2 mm	S-BMCA8020
Mini Couple Abutment for Short Implant	ø 8.0	3 mm	S-BMCA8030

#### ANGLED ABUTMENTS

Mini Angled Abutment	ø 3.7	1 mm 15°	S-BMAA37151
Mini Angled Abutment	ø 3.7	2 mm 15°	S-BMAA37152
Mini Angled Abutment	ø 3.7	4 mm 15°	S-BMAA37154

Angled Abutment	ø 4.5	2 mm 15°	S-BAA45152
Angled Abutment	ø 4.5	2 mm 25°	S-BAA45252
Angled Abutment	ø 4.5	4 mm 15°	S-BAA45154
Angled Abutment	ø 4.5	4 mm 25°	S-BAA45254

Angled Abutment	ø 5.5	2 mm 15°	S-BAA55152
Angled Abutment	ø 5.5	2 mm 25°	S-BAA55252
Angled Abutment	ø 5.5	4 mm 15°	S-BAA55154
Angled Abutment	ø 5.5	4 mm 25°	S-BAA55254

Angled Abutment	ø 6.5	2 mm 15°	S-BAA65152
Angled Abutment	ø 6.5	2 mm 25°	S-BAA65252
Angled Abutment	ø 6.5	4 mm 15°	S-BAA65154
Angled Abutment	ø 6.5	4 mm 25°	S-BAA65254

#### TEMPORARY ABUTMENTS

Mini Temporary Abutment	ø 3.7		S-BMTA37
Temporary Abutment	ø 4.5		S-BTA45
Temporary Abutment	ø 5.5		S-BTA55

#### IMPRESSION COPING PICK UP

Mini Impression Coping Pick up	ø 3.7		S-BMICP37L
Mini Impression Coping Pick up	ø 3.7		S-BMICP37S

Impression Coping Pick up	ø 4.5		S-BICP45L
Impression Coping Pick up	ø 4.5		S-BICP45S

Impression Coping Pick up	ø 5.5		S-BICP55L
Impression Coping Pick up	ø 5.5		S-BICP55S

Impression Coping Pick up	ø 6.5		S-BICP65L
Impression Coping Pick up	ø 6.5		S-BICP65S

Impression Coping Pick up for short	ø 5.0		S-BMSICP50
Impression Coping Pick up for short	ø 6.0		S-BMSICP60

PRODUCT	PLATFORM	HEIGHT	CODE
<b>IMPRESSION COPING TRANSFER</b>			
Mini Impression Coping Transfer	ø 3.7		S-BMICT37L
Mini Impression Coping Transfer	ø 3.7		S-BMICT37S
Mini Impression Coping Transfer for short	ø 5.0		S-BMSICT50
Mini Impression Coping Transfer for short	ø 6.0		S-BMSICT60
Impression Coping Transfer	ø 4.5		S-BIP45

#### FIXTURE LAB ANALOG

Mini Fixture Lab Analog	ø 3.3	12 mm	S-BMFLA
Fixture Lab Analog	ø 4.3	12 mm	S-BFLA

#### BONE LEVEL SOLID SYSTEM

##### SOLID ABUTMENTS

Mini Solid Abutment	ø 3.7	1 mm	S-BMSA3710
Mini Solid Abutment	ø 3.7	2 mm	S-BMSA3720
Mini Solid Abutment	ø 3.7	3 mm	S-BMSA3730
Mini Solid Abutment	ø 3.7	4 mm	S-BMSA3740
Mini Solid Abutment	ø 3.7	5 mm	S-BMSA3750

Solid Abutment	ø 4.5	1 mm	S-BSA4510
Solid Abutment	ø 4.5	2 mm	S-BSA4520
Solid Abutment	ø 4.5	3 mm	S-BSA4530
Solid Abutment	ø 4.5	4 mm	S-BSA4540
Solid Abutment	ø 4.5	5 mm	S-BSA4550

Solid Abutment	ø 5.5	1 mm	S-BSA5510
Solid Abutment	ø 5.5	2 mm	S-BSA5520
Solid Abutment	ø 5.5	3 mm	S-BSA5530
Solid Abutment	ø 5.5	4 mm	S-BSA5540
Solid Abutment	ø 5.5	5 mm	S-BSA5550

#### SOLID ABUTMENT IMPRESSION CAPS

Mini Solid Abutment Impression Cap	ø 3.7		S-BMSAICap37
Solid Abutment Impression Cap	ø 4.5		S-BSAICap45
Solid Abutment Impression Cap	ø 5.5		S-BSAICap55

#### SOLID LAB ANALOG

Mini Solid Lab Analog	ø 3.7		S-BMSALA37
Solid Lab Analog	ø 4.5		S-BSALA45
Solid Lab Analog	ø 5.5		S-BSALA55

#### BONE LEVEL OCTA SYSTEM

##### OCTA ABUTMENTS

Mini Octa Abutment	ø 4.8	1 mm	S-BMOA10
Mini Octa Abutment	ø 4.8	2 mm	S-BMOA20
Mini Octa Abutment	ø 4.8	3 mm	S-BMOA30
Mini Octa Abutment	ø 4.8	4 mm	S-BMOA40
Mini Octa Abutment	ø 4.8	5 mm	S-BMOA50

PRODUCT	PLATFORM	HEIGHT	CODE
Octa Abutment	∅ 4.8	1 mm	S-BOA10
Octa Abutment	∅ 4.8	2 mm	S-BOA20
Octa Abutment	∅ 4.8	3 mm	S-BOA30
Octa Abutment	∅ 4.8	4 mm	S-BOA40
Octa Abutment	∅ 4.8	5 mm	S-BOA50

#### OCTA ABUTMENT HEALING CAP

Octa Abutment Healing Cap	∅ 4.9	4 mm	S-OAHCap
Octa Abutment Healing Cap Large	∅ 4.9	6 mm	S-OAHCAP-L

#### OCTA SYSTEM PLASTIC CYLINDER OCTA

Octa System Plastic Cylinder Octa	∅ 5.2 mm		S-OAPICy
-----------------------------------	----------	--	----------

#### OCTA SYSTEM PLASTIC CYLINDER NON-OCTA

Octa System Plastic Cylinder Non-Octa	∅ 5.2 mm		S-OAPICyN
---------------------------------------	----------	--	-----------

#### OCTA SYSTEM TITANIUM CYLINDER

Octa System Plastic Cylinder	∅ 5.2 mm		S-OATICy
------------------------------	----------	--	----------

#### OCTA SYSTEM TITANIUM CYLINDER NON-OCTA

Octa System Plastic Cylinder Non-Hex	∅ 5.2		S-OATICyN
--------------------------------------	-------	--	-----------

#### OCTA SYSTEM IMPRESSION COPING PICK UP

Octa System Impression Coping Pick Up	∅ 5.5	10 mm	S-ICOAP
---------------------------------------	-------	-------	---------

#### OCTA SYSTEM IMPRESSION COPING PICK UP NON-OCTA

Octa System Impression Coping Pick Up Non-Hex	∅ 5.5	10 mm	S-ICOAPN
---	-------	-------	----------

#### OCTA SYSTEM IMPRESSION COPING TRANSFER

Octa System Impression Coping Transfer	∅ 5.0	7 mm	S-ICOAT
--	-------	------	---------

#### OCTA SYSTEM IMPRESSION COPING TRANSFER NON-OCTA

Octa System Impression Coping Transfer Non-Hex	∅ 5.0	7 mm	S-ICOATN
--	-------	------	----------

#### OCTA ABUTMENT LAB ANALOG

Octa Abutment Lab Analog	∅ 4.8	14 mm	S-OALA
--------------------------	-------	-------	--------

### BONE LEVEL O-RING SYSTEM

#### O-RING ABUTMENTS

Mini O-Ring Abutment	∅ 2.9	0.5 mm	S-BMORA05
Mini O-Ring Abutment	∅ 3.7	2 mm	S-BMORA20
Mini O-Ring Abutment	∅ 3.7	4 mm	S-BMORA40

O-Ring Abutment	∅ 3.42	0.5 mm	S-BORA05
O-Ring Abutment	∅ 3.42	1 mm	S-BORA10

O-Ring Abutment	∅ 4.5	2 mm	S-BORA20
O-Ring Abutment	∅ 4.5	3 mm	S-BORA30
O-Ring Abutment	∅ 4.5	4 mm	S-BORA40
O-Ring Abutment	∅ 4.5	6 mm	S-BORA60

PRODUCT	PLATFORM	HEIGHT	CODE
<b>DIRECTIONAL RINGS</b>			
Directional Rings			S-ORRa

#### RETAINERS

Retainer	∅ 5.5	2 mm	S-ORR
Retainer	∅ 5.5	3.5 mm	S-ORR-C
Retainer Pink Cap			S-ORG-PCap
Retainer Yellow Cap			S-ORG-YCap
Retainer Green Cap			S-ORG-GCap

#### O-RING LAB ANALOG

O-Ring Lab Analog	∅ 4.1		S-ORLA
-------------------	-------	--	--------

### BONE LEVEL LOCATOR SYSTEM

#### LOCATOR ABUTMENTS

Locator Mini Abutment	∅ 3.86mm	1.5mm	S-BMLA15
Locator Mini Abutment	∅ 3.86mm	2.5mm	S-BMLA25
Locator Mini Abutment	∅ 3.86mm	3.5mm	S-BMLA35
Locator Mini Abutment	∅ 3.86mm	4.5mm	S-BMLA45
Locator Mini Abutment	∅ 3.86mm	5.5mm	S-BMLA55
Locator Mini Abutment	∅ 3.86mm	6.5mm	S-BMLA65

Locator Abutment	∅ 3.86mm	1.5mm	S-BLA15
Locator Abutment	∅ 3.86mm	2.5mm	S-BLA25
Locator Abutment	∅ 3.86mm	3.5mm	S-BLA35
Locator Abutment	∅ 3.86mm	4.5mm	S-BLA45
Locator Abutment	∅ 3.86mm	5.5mm	S-BLA55
Locator Abutment	∅ 3.86mm	6.5mm	S-BLA65

#### LOCATOR MATRIX

Locator Matrix (includes processing cap)	∅ 5.5mm		S-LMHt
Locator Cap (includes 4 pcs caps)	∅ 4.86 mm		S-LRI

#### LOCATOR ANALOG

Locator Lab Analog	∅ 4 mm		S-BLLA
--------------------	--------	--	--------

### BONE LEVEL MULTI-UNIT SYSTEM

#### MULTI-UNIT ANGLED ABUTMENT

Bone Level Multi Unit Angled Abutment	∅ 4.8x17	1.0 mm	S-BMUA48171
Bone Level Multi Unit Angled Abutment	∅ 4.8x17	2.5 mm	S-BMUA48172
Bone Level Multi Unit Angled Abutment	∅ 4.8x17	3.5 mm	S-BMUA48173
Bone Level Multi Unit Angled Abutment	∅ 4.8x35	1.0 mm	S-BMUA48351
Bone Level Multi Unit Angled Abutment	∅ 4.8x35	2.5 mm	S-BMUA48352
Bone Level Multi Unit Angled Abutment	∅ 4.8x35	3.5 mm	S-BMUA48353

#### MULTI-UNIT ABUTMENT

Bone Level Multi Unit Mini Abutment	∅ 4.8	0.5 mm	S-BMMUA4805
Bone Level Multi Unit Mini Abutment	∅ 4.8	1.5 mm	S-BMMUA4815
Bone Level Multi Unit Mini Abutment	∅ 4.8	2.5 mm	S-BMMUA4825

PRODUCT	PLATFORM	HEIGHT	CODE
Bone Level Multi Unit Mini Abutment	ø 4.8	3.5 mm	S-BMMUA4835
Bone Level Multi Unit Mini Abutment	ø 4.8	4.5 mm	S-BMMUA4845
Bone Level Multi Unit Abutment	ø 4.8	0.5 mm	S-BMUA4805
Bone Level Multi Unit Abutment	ø 4.8	1.5 mm	S-BMUA4815
Bone Level Multi Unit Abutment	ø 4.8	2.5 mm	S-BMUA4825
Bone Level Multi Unit Abutment	ø 4.8	3.5 mm	S-BMUA4835
Bone Level Multi Unit Abutment	ø 4.8	4.5 mm	S-BMUA4845
<b>MULTI UNIT HEALING CAP</b>			
Multi Unit Healing Cap	ø 4.8	4.0 mm	S-MUAHCAP
Multi Unit Healing Cap	ø 4.8	6.0 mm	S-MUAHCAPL
<b>MULTI UNIT IMPRESSION COPING PICK UP</b>			
Bone Level Multi Unit Impression Coping Pick Up	ø 4.8		S-BMUAICP
<b>MULTI UNIT ABUTMENT PLASTIC CYLINDER</b>			
Multi Unit Abutment Plastic Cylinder	ø 4.8		S-MUAPICyN
<b>MULTI UNIT ABUTMENT TITANIUM CYLINDER</b>			
Multi Unit Abutment Titanium Cylinder	ø 4.8		S-MUATICyN
<b>MULTI UNIT LAB ANALOG</b>			
Multi Unit Abutment Lab Analog	ø 4.8	14 mm	S-MUALA
<b>BAR SYSTEM</b>			
Bar burnout			S-IB
Bar Positioning			S-IBP
Bar Housing			S-IBH
Bar Clip Yellow			S-IBClip-y
Bar Clip Pink			S-IBClip-P
Bar Clip Red			S-IBClip-r
Tool for Inserting Clip			S-BMIA
Key for Parallelometer			S-TFIC
<b>TISSUE LEVEL FIXTURES</b>			
<b>TISSUE LEVEL FIXTURE HYBRID</b>			
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	8 mm	S-TFHR3708
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	10 mm	S-TFHR3710
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	12 mm	S-TFHR3712
Tissue Level Wide Fixture Hybrid SRA	ø 3.7	14 mm	S-TFHR3714
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	8 mm	S-TFHR4308
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	10 mm	S-TFHR4310
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	12 mm	S-TFHR4312
Tissue Level Wide Fixture Hybrid SRA	ø 4.3	14 mm	S-TFHR4314
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	8 mm	S-TFHR4808
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	10 mm	S-TFHR4810

PRODUCT	PLATFORM	HEIGHT	CODE
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	12 mm	S-TFHR4812
Tissue Level Wide Fixture Hybrid SRA	ø 4.8	14 mm	S-TFHR4814
<b>TISSUE LEVEL SHORT FIXTURE STRAIGHT</b>			
Tissue Level Short Fixture Straight SRA	ø 4.3	4 mm	S-TFSR4304
Tissue Level Short Fixture Straight SRA	ø 4.3	5 mm	S-TFSR4305
Tissue Level Short Fixture Straight SRA	ø 4.3	6 mm	S-TFSR4306
Tissue Level Short Fixture Straight SRA	ø 4.8	4 mm	S-TFSR4804
Tissue Level Short Fixture Straight SRA	ø 4.8	5 mm	S-TFSR4805
Tissue Level Short Fixture Straight SRA	ø 4.8	6 mm	S-TFSR4806
Tissue Level Short Fixture Straight SRA	ø 5.5	4 mm	S-TWFSR5504
Tissue Level Short Fixture Straight SRA	ø 5.5	5 mm	S-TWFSR5505
Tissue Level Short Fixture Straight SRA	ø 5.5	6 mm	S-TWFSR5506
<b>TISSUE LEVEL WIDE FIXTURE HYBRID</b>			
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	8 mm	S-TWFHR5508
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	10 mm	S-TWFHR5510
Tissue Level Wide Fixture Hybrid SRA	ø 5.5	12 mm	S-TWFHR5512
<b>TISSUE LEVEL HEALING ABUTMENTS</b>			
Tissue Level Healing Abutment	ø 5.2	3.5 mm	S-THA35
Tissue Level Healing Abutment	ø 5.2	5.5 mm	S-THA55
Tissue Level Healing Abutment	ø 6.4	3.5 mm	S-THAW35
Tissue Level Healing Abutment	ø 6.4	5.5 mm	S-THAW55
<b>TISSUE LEVEL COUPLE SYSTEM</b>			
<b>TISSUE LEVEL COUPLE ABUTMENTS</b>			
Tissue Level Couple Abutment	ø 5.2	1 mm	S-TCA10
Tissue Level Couple Abutment	ø 5.2	2 mm	S-TCA20
Tissue Level Couple Abutment	ø 5.2	3 mm	S-TCA30
Tissue Level Couple Abutment	ø 5.2	4 mm	S-TCA40
Tissue Level Couple Abutment Wide	ø 6.4	1 mm	S-TCAW10
Tissue Level Couple Abutment Wide	ø 6.4	2 mm	S-TCAW20
Tissue Level Couple Abutment Wide	ø 6.4	3 mm	S-TCAW30
Tissue Level Couple Abutment Wide	ø 6.4	4 mm	S-TCAW40
<b>TISSUE LEVEL ANGLED ABUTMENTS</b>			
Tissue Level Angled Abutment	ø 3.7	15°	S-TAA15
Tissue Level Angled Abutment	ø 3.7	25°	S-TAA25
<b>TISSUE LEVEL SYNOCTA ABUTMENT OCTA</b>			
Tissue Level Synocta Abutment Octa	ø 3.5	6 mm	S-TSOA
<b>TISSUE LEVEL IMPRESSION COPING PICK UP OCTA</b>			
Tissue Level Impression Coping Pick Up Octa	ø 5.5	11.75 mm	S-TICP

PRODUCT	PLATFORM	HEIGHT	CODE
Tissue Level Impression Coping Pick Up Wide Octa	ø 6.4	11.75 mm	S-TICPW
<b>TISSUE LEVEL IMPRESSION COPING TRANSFER</b>			
Tissue Level Impression Coping Transfer Octa	ø 5.5	9 mm	S-TICT
<b>TISSUE LEVEL FIXTURE LAB ANALOG</b>			
Tissue Level Fixture Lab Analog	ø 4.8	14 mm	S-TFLA
Tissue Level Fixture Large Lab Analog	ø 6.0	14 mm	S-TFLLA
<b>TISSUE LEVEL SOLID SYSTEM</b>			
<b>TISSUE LEVEL SOLID ABUTMENTS</b>			
Tissue Level Solid Abutment	ø 3.5	4 mm	S-TSA40
Tissue Level Solid Abutment	ø 3.5	5.5 mm	S-TSA55
Tissue Level Solid Abutment	ø 3.5	7 mm	S-TSA70
Tissue Level Solid Abutment Large	ø 4.5	4 mm	S-TSAL40
Tissue Level Solid Abutment Large	ø 4.5	5 mm	S-TSAL50
Tissue Level Solid Abutment Large	ø 4.5	7 mm	S-TSAL70
<b>TISSUE LEVEL SOLID ABUTMENT POSITIONING CYLINDER</b>			
Tissue Level Solid Abutment Positioning Cylinder	ø 10.2		S-TSAPoCy
<b>TISSUE LEVEL SOLID ABUTMENT PLASTIC COPING</b>			
Tissue Level Solid Plastic Coping Crown	ø 5.3	10 mm	S-TSAPICc
Tissue Level Solid Plastic Coping Bridge	ø 5.3	10 mm	S-TSAPICb
<b>TISSUE LEVEL SOLID ABUTMENT IMPRESSION CAP</b>			
Tissue Level Solid Abutment Impression Cap	ø 10.2		S-TSAICap
<b>TISSUE LEVEL SOLID ABUTMENT LAB ANALOG</b>			
Tissue Level Solid Abutment Lab Analog	ø 4.8	14 mm	S-TSALA40
Tissue Level Solid Abutment Lab Analog	ø 4.8	15.5 mm	S-TSALA55
Tissue Level Solid Abutment Lab Analog	ø 4.8	17 mm	S-TSALA70
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	14 mm	S-TSALLA40
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	15.5 mm	S-TSALLA55
Tissue Level Solid Abutment Large Lab Analog	ø 6.0	17 mm	S-TSALLA70
<b>TISSUE LEVEL OCTA SYSTEM</b>			
<b>TISSUE LEVEL OCTA ABUTMENT</b>			
Tissue Level Octa Abutment	ø 3.5		S-TOA
<b>OCTA ABUTMENT HEALING CAP</b>			
Octa Abutment Healing Cap	ø 4.9	4 mm	S-OAHCAP
Octa Abutment Healing Cap Large	ø 4.9	6 mm	S-OAHCAP-L
<b>TISSUE LEVEL OCTA SYSTEM PLASTIC CYLINDER</b>			
Octa System Plastic Cylinder Octa	ø 5		S-OAPICy

PRODUCT	PLATFORM	HEIGHT	CODE
<b>TISSUE LEVEL OCTA ABUTMENT PLASTIC CYLINDER NON-OCTA</b>			
Tissue Level Octa Abutment Plastic Cylinder Non-Octa	ø 5		S-OAPICyN
<b>TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING PICK UP</b>			
Tissue Level Octa System Impression Coping Pick Up Octa	ø 5.5		S-ICOAP
Tissue Level Octa System Impression Coping Pick Up Non-Octa	ø 5.5		S-ICOAPN
<b>TISSUE LEVEL OCTA SYSTEM IMPRESSION COPING TRANSFER</b>			
Tissue Level Octa System Impression Coping Transfer Octa	ø 5.0		S-ICOAT
Tissue Level Octa System Impression Coping Transfer Non-Octa	ø 5.0		S-ICOATN
<b>TISSUE LEVEL OCTA ABUTMENT LAB ANALOG</b>			
Tissue Level Octa Abutment Lab Analog	ø 4.8		S-OALA
<b>TISSUE LEVEL O-RING SYSTEM</b>			
<b>TISSUE LEVEL O-RING ABUTMENTS</b>			
Tissue Level O-Ring Abutment	ø 3.5	0.5 mm	S-TORA05
Tissue Level O-Ring Abutment	ø 3.5	2 mm	S-TORA20
Tissue Level O-Ring Abutment	ø 3.5	4 mm	S-TORA40
<b>DIRECTIONAL RINGS</b>			
Directional Rings			S-ORRa
<b>RETAINERS</b>			
Retainer	ø 5.5	2 mm	S-ORR
Retainer	ø 5.5	3.5 mm	S-ORR-C
Retainer Pink Cap			S-ORG-PCap
Retainer Yellow Cap			S-ORG-YCap
Retainer Green Cap			S-ORG-GCap
<b>TISSUE LEVEL O-RING LAB ANALOG</b>			
Tissue Level O-Ring Lab Analog	ø 4.1		S-ORLA
<b>TISSUE LEVEL LOCATOR SYSTEM</b>			
<b>TISSUE LEVEL LOCATOR ABUTMENTS</b>			
Tissue Level Locator Abutment	ø 4.9	1.5 mm	S-TLA15
Tissue Level Locator Abutment	ø 4.9	2.5 mm	S-TLA25
Tissue Level Locator Abutment	ø 4.9	3.5 mm	S-TLA35
Tissue Level Locator Abutment	ø 4.9	4.5 mm	S-TLA45
Tissue Level Locator Abutment	ø 4.9	5.5 mm	S-TLA55
Tissue Level Locator Abutment	ø 4.9	6.5 mm	S-TLA65
<b>LOCATOR MATRIX</b>			
Locator Matrix (includes processing cap)	ø 5.5mm		S-LMHt
Locator Cap (includes 4 pcs caps)	ø 4.86 mm		S-LRI
<b>TISSUE LEVEL LOCATOR LAB ANALOG</b>			
Tissue Level Locator Lab Analog			S-TLLA

PRODUCT	PLATFORM	HEIGHT	CODE
---------	----------	--------	------

### ZYGOMA IMPLANT

#### ZYGOMA FIXTURE

Zygoma Fixture	ø 4.3	35 mm	S-ZYG35
Zygoma Fixture	ø 4.3	37.5 mm	S-ZYG37.5
Zygoma Fixture	ø 4.3	40 mm	S-ZYG40
Zygoma Fixture	ø 4.3	42.5 mm	S-ZYG42.5
Zygoma Fixture	ø 4.3	45 mm	S-ZYG45
Zygoma Fixture	ø 4.3	47.5 mm	S-ZYG47.5
Zygoma Fixture	ø 4.3	50 mm	S-ZYG50
Zygoma Fixture	ø 4.3	52.5 mm	S-ZYG52.5
Zygoma Fixture	ø 4.3	55 mm	S-ZYG55
Zygoma Fixture	ø 4.3	57.5 mm	S-ZYG57.5
Zygoma Fixture	ø 4.3	60 mm	S-ZYG60

#### ATROPHIC MAXILLARY ZYGOMA FIXTURE

Atrophic Maxillary Zygoma Fixture	ø 4.3	35 mm	S-ZAM35
Atrophic Maxillary Zygoma Fixture	ø 4.3	37.5 mm	S-ZAM37.5
Atrophic Maxillary Zygoma Fixture	ø 4.3	40 mm	S-ZAM40
Atrophic Maxillary Zygoma Fixture	ø 4.3	42.5 mm	S-ZAM42.5
Atrophic Maxillary Zygoma Fixture	ø 4.3	45 mm	S-ZAM45
Atrophic Maxillary Zygoma Fixture	ø 4.3	47.5 mm	S-ZAM47.5
Atrophic Maxillary Zygoma Fixture	ø 4.3	50 mm	S-ZAM50
Atrophic Maxillary Zygoma Fixture	ø 4.3	52.5 mm	S-ZAM52.5
Atrophic Maxillary Zygoma Fixture	ø 4.3	55 mm	S-ZAM55
Atrophic Maxillary Zygoma Fixture	ø 4.3	57.5 mm	S-ZAM57.5
Atrophic Maxillary Zygoma Fixture	ø 4.3	60 mm	S-ZAM60

#### ONCOLOGICAL ZYGOMA FIXTURE

Oncological Zygoma Fixture	ø 4.3	25 mm	SZ-ONC25
Oncological Maxillary Zygoma Fixture	ø 4.3	27.5 mm	SZ-ONC27.5
Oncological Maxillary Zygoma Fixture	ø 4.3	30 mm	SZ-ONC30
Oncological Maxillary Zygoma Fixture	ø 4.3	32.5 mm	SZ-ONC32.5

#### ABUTMENTS

Zygoma Abutment		3 mm	S-ZCA4330
Zygoma Abutment		5 mm	S-ZCA4350

#### ZYGOMA ANGLED ABUTMENTS

Zygoma Angled Abutment		3 mm	S-ZAA43203
Zygoma Angled Abutment		5 mm	S-ZAA43205

#### ZYGOMA MULTI UNIT ABUTMENT

Zygoma Multi Unit Abutment		0.5 mm	S-ZMUA431705
Zygoma Multi Unit Abutment		1 mm	S-ZMUA43171
Zygoma Multi Unit Abutment		2 mm	S-ZMUA43172
Zygoma Multi Unit Abutment		3 mm	S-ZMUA43173
Zygoma Multi Unit Abutment		5 mm	S-ZMUA43175

PRODUCT	PLATFORM	HEIGHT	CODE
---------	----------	--------	------

### CAPS

Zygoma Healing Cap		3 mm	S-ZHA4330
Zygoma Healing Cap		5 mm	S-ZHA4350

### COPING

Zygoma Impression Coping Pick Up			S-ZICP
----------------------------------	--	--	--------

### PLASTIC CYLINDER

Zygoma Plastic Cylinder			S-ZPICYN
-------------------------	--	--	----------

### LAB ANALOG

Zygoma Lab Analog			S-ZLA
-------------------	--	--	-------

### DRIVERS

Implant Drivers (For ratchet and contra angle)

Bone Level Mini Fixture Driver Machine Short			S-BMFDMS
Bone Level Mini Fixture Driver Machine Long			S-BMFDML

Bone Level Mini Fixture Driver Ratchet Short			S-BMFDMS
Bone Level Mini Fixture Driver Ratchet Long			S-BMFDRL

Bone Level Fixture Driver Machine Short			S-BFDMS
Bone Level Fixture Driver Machine Long			S-BFDML

Bone Level Fixture Driver Ratchet Short			S-BFDMS
Bone Level Fixture Driver Ratchet Long			S-BFDRL

Tissue Level Fixture Driver Ratchet Short			S-TFDMS
Tissue Level Fixture Driver Ratchet Long			S-TFDRL

Tissue Level Fixture Driver Machine Short			S-TFDMS
Tissue Level Fixture Driver Machine Long			S-TFDML

Hex driver (extra short, short, long, extra long)			S-HDRS / S-HDRL
---	--	--	-----------------

Hex Driver Contra Angle (extra short, short, long, extra long)			S-HDRS / S-HDRL
--	--	--	-----------------

O-ring Driver			S-ORADS/L
---------------	--	--	-----------

Octa Driver			S-OADS/L
-------------	--	--	----------

Solid Driver			S-SAD
--------------	--	--	-------

ORRO Insert Adaptor			S-OIA
---------------------	--	--	-------

Bar Matrix Insert Adaptor			S-BMIA
---------------------------	--	--	--------



PRODUCT	PLATFORM	HEIGHT	CODE
<b>SCREWS</b>			
<b>CLOSING SCREW</b>			
<b>BONE LEVEL</b>			
Bone Level Mini Closing Screw	∅ 2.8	4.5 mm	S-BMCS
Bone Level Closing Screw	∅ 3.36	6.5 mm	S-BCS
<b>TISSUE LEVEL</b>			
Tissue Level Closing-Screw Small	∅ 3.5	5.35 mm	S-TCSS
Tissue Level Closing-Screw Large	∅ 5.2	6.4mm	S-TCSL
Tissue Level Closing-Screw Wide	∅ 6.4	6.4 mm	S-TCSW
<b>ABUTMENT SCREWS</b>			
<b>BONE LEVEL</b>			
Bone Level Mini Abutment Screw	∅ 2.1	7.64 mm	S-BMAS
Bone Level Abutment Screw	∅ 2.3	9.32 mm	S-BAS
Bone Level Multi-Unit Abutment Screw	∅ 1.96	9 mm	S-BMUAS
<b>TISSUE LEVEL</b>			
Tissue Level Abutment Screw	∅ 2.5	7.1	S-TAS
<b>COMMON SCREWS</b>			
Octa Abutment Screw	∅ 2.5	4.5 mm	S-OAS
Multi-Unit Abutment Healing Screw	∅ 2.2	4 mm	S-MUAHS
<b>DIGITAL SOLUTIONS</b>			
<b>BONE LEVEL SCAN BODY</b>			
Bone Level Mini Scan Body			S-PMSB
Bone Level Scan Body			S-PSB
<b>BONE LEVEL DIGITAL ANALOG</b>			
Bone Level Mini Digital Analog			S-BMFDLA
Bone Level Digital Analog			BFDLA
<b>BONE LEVEL PREMILL</b>			
Bone Level Mini Premill	10 mm		S-BMPA10
Bone Level Mini Premill	14 mm		S-BMPA14
Bone Level Premill	10 mm		S-BPA10
Bone Level Premill	14 mm		S-BPA14
Bone Level Premill	16 mm		S-BPA16
<b>TISSUE LEVEL SCAN BODY</b>			
Tissue Level Mini Scan Body			S-TSB
Tissue Level Scan Body			S-TWFDLA

PRODUCT	PLATFORM	HEIGHT	CODE
<b>TISSUE LEVEL DIGITAL ANALOG</b>			
Tissue Level Mini Digital Analog			S-TFDLA
Tissue Level Digital Analog			S-TWFDLA
<b>TISSUE LEVEL PREMILL</b>			
Tissue Level Mini Premill	10 mm		S-TPA10
Tissue Level Mini Premill	14 mm		S-TPA14
Tissue Level Mini Premill	16 mm		S-TPA16
Tissue Level Premill	10 mm		S-TWPA10
Tissue Level Premill	14 mm		S-TWPA14
Tissue Level Premill	16 mm		S-TWPA16
<b>MULTI UNIT ABUTMENT</b>			
Multi Unit Abutment Scan Body			S-PMUSB
Multi Unit Abutment Digital Analog			S-BMDLA
Multi Unit Abutment ZR Base			S-MUTB
<b>OCTA ABUTMENT</b>			
Octa Abutment Scan Body			S-POSB
Octa Abutment Digital Analog			S-BODLA
Octa Abutment ZR Base			S-OTB
<b>Ti-BASE ABUTMENT</b>			
Ti-Base Mini Abutment			S-BMTIA05
Ti-Base Abutment			S-BTIA05
<b>Ti-BASE SCAN BODY</b>			
Ti-Base Mini Scan Body			S-TIMSB
Ti-Base Scan Body			S-TISB



