


» Multi Unit Abutment

ABUTMENT OPTIONS FOR ALL ON FOUR TECHNIQUE



ALL ON FOUR TECHNIQUE

It is a dental implant application technique preferred for fully edentulous patients. The dental prosthesis is fixed by placing 4 implants at certain angles to the patient, thus providing a dental solution at the most appropriate cost and time without the need for an additional operation.

implant**swiss**
IMPLANT TECHNOLOGY
SMART
CHOICE

MULTI UNIT ABUTMENTS



ø4.8

only for ø3.3 implant

0.5 mm S-BMMUA4805
1.5 mm S-BMMUA4815
2.5 mm S-BMMUA4825
3.5 mm S-BMMUA4835
4.5 mm S-BMMUA4845



ø4.8

0.5 mm S-BMUA4805
1.5 mm S-BMUA4815
2.5 mm S-BMUA4825
3.5 mm S-BMUA4835
4.5 mm S-BMUA4845

System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » : Abutment Screw
25 Ncm Oclusal Screw 20 Ncm
Driver » 1,2 Hex Driver

MULTI UNIT ANGLED ABUTMENTS



ø4.8 (17°)

1.0 mm S-BMUA48171
2.5 mm S-BMUA48172
3.5 mm S-BMUA48173



ø4.8 (35°)

1.0 mm S-BMUA48351
2.5 mm S-BMUA48352
3.5 mm S-BMUA48353

System » Bone Level
Material » Ti Grade 5 ELI
Torque Value » : Abutment Screw
25 Ncm Oclusal Screw 20 Ncm
Driver » 1,2 Hex Driver

MULTI UNIT ABUTMENT HEALING CAP



ø4.8

S-BMUAHCap

System » Multi Unit System
Material » Ti Grade 5 ELI

MULTI UNIT IMPRESSION COPING PICK UP



ø4.8

S-BMUAICP

System » Multi Unit System
Material » Ti Grade 5 ELI

MULTI UNIT ABUTMENT PLASTIC CYLINDER



S-MUAPICY

System » Multi Unit System
Material » POM

MULTI UNIT ABUTMENT TITANIUM CYLINDER



S-MUATICy

System » Multi Unit System
Material » Ti Grade 5 ELI

MULTI UNIT LAB ANALOG



S-MUALA

ø 4.8×14mm

System » Bone Level
Material » Ti Grade 5 ELI